

## Toluene, Blood

Order Name: **TOLUENE**  
Test Number: 3610275  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Toluene, Blood	Headspace Gas Chromatography/Flame Ionization Detection (GC-FID)	9637-0

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	20 mL (7 mL)	Whole Blood	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated
Instructions	<p><b>Specimen Type:</b> Green-top (heparin) tube  <b>Specimen Storage:</b> Refrigerate; specimen may also be frozen.  <b>Specimen Collection:</b> Collect adequate number of tubes to obtain the required specimen volume. Tubes should be completely filled to prevent loss of volatile compound into headspace. Collect prior to last shift of workweek.  <b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc TC 70123</p>			

GENERAL INFORMATION	
Expected TAT	10 - 14 days
Performing Labcorp Test Code	790580
Notes	Labcorp Test Code: 790580
CPT Code(s)	84600
Lab Section	Reference Lab