

Fanconi Anemia, Group C (FANCC) 2 Mutations

Order Name: **Fanconi Anem Mut**

Test Number: 6905009

Revision Date: 04/15/2017

TEST NAME	METHODOLOGY	LOINC CODE
Fanconi Anemia, Group C (FANCC) 2 Mutations	Polymerase Chain Reaction	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	7 mL (3 mL)	Whole Blood	EDTA (Lavender Top)	Room Temperature
Instructions	<p>Notes: 3 mL whole blood, 5 mL amniotic fluid, 10 mg CVS, OR two buccal swabs</p> <p>Specimen Type: Lavender-top (EDTA) tube, yellow-top (ACD) tube, sterile plastic conical tube or two confluent T-25 flasks for fetal testing, OR LabCorp buccal swab kit</p> <p>Specimen Storage: Maintain specimen at room temperature or refrigerate.</p> <p>Specimen Collection: Not Available</p> <p>Special Instructions: If cultured cells are needed, an additional 7-12 days may be required. Additional culture fee may be included.</p>			

GENERAL INFORMATION	
Expected TAT	8-15 days
Performing Labcorp Test Code	511212
Notes	Labcorp Test Code: 511212
CPT Code(s)	81242
Lab Section	Reference Lab