

## ENT Allergy Panel 1 (Large panel)

Order Name: **ENT PNL 1**  
Test Number: 6908087  
Revision Date: 08/01/2023

TEST NAME	METHODOLOGY	LOINC CODE
ENT Allergy Panel 1 (Large panel)	ImmunoCAP	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	18 mL (9.4 mL)	Serum	Clot Activator (SST or Red No-Gel)	Room Temperature
Instructions	<p><b>Specimen Type:</b> Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container.</p> <p><b>Stability Requirements:</b> Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)</p>			

### GENERAL INFORMATION

Expected TAT	3-5 days
--------------	----------

Notes			
	Alternaria alternata Allergen	Cottonwood Tree Allergen	Mugwort Allergen
	Aspergillus fumigatus Allergen	Dandelion Allergen	Mulberry Tree White Allergen
	Bahia Grass Allergen	Dog Dander Allergen	Nettle Allergen
	Bermuda Grass Allergen	Dust Mite (D. farinae) Allergen	Oak Tree Allergen
	Bipolaris spicifera IgE	Elm Tree Allergen	Pecan, Hickory Tree Allergen
	Box Elder (Maple) Tree Allergen	Epicoccum purpurascens Allergen	Penicillin Chrysogenium Allergen
	Cat Dander Allergen	Feathers Allergen Panel	Plantain (English), Ribwort Allergen
	Cedar/Red Tree IgE	Firebush (Kochia) Allergen	Privet Tree Allergen
	Cephalosporium IgE	Giant Ragweed Allergen	Ragweed Common Short Allergen
	Chaetomium globosum IgE	Goosefoot, Lamb's Quarters Allergen	Rough Marshelder Allergen
	Cladosporium herbarum Allergen	House Dust Mite Allergen	Saltwort, Russian Thistle Allergen
	Cocklebur Allergen	Immunoglobulin IgE	Sheep Sorrel Allergen
	Cockroach, American IgE	Johnson Grass Allergen	Timothy Grass Allergen
	Cockroach, German Allergen	Mesquite Tree Allergen	White Ash Tree Allergen
	Common Pigweed Allergen	Mountain Juniper Tree Allergen	White Pine Tree IgE
	Common Silver Birch Tree Allergen	Mucor racemosus Allergen	

Allergens performed at Labcorp Phoenix

CPT Code(s)	86003 x47
Lab Section	Reference Lab

